

H100B

Pulse Oximeter



EDAN

A world of potential

H100B Pulse Oximeter

As a handheld pulse oximeter, H100B is designed to serve the needs of both long-term monitoring and spot check. It's dedicated to offering reliable, accurate and sensitive measurement of blood oxygen saturation under even harsh conditions.

Measured Parameters: SpO₂ + PR

- LCD display with backlight control
- Dual work modes: monitoring & spot check
- Numeric display with plethysmogram display
- Trend review
- Pulse-tone modulation (Pitch Tone)
- Powerful data storage capacity (up to 300 hours)
- PatientCare Viewer software for PC data management
- 4 x AA batteries for up to 48 hours of work

300 h

SpO₂ Review

100

Patient ID

- Rechargeable battery for up to 30 hours of work
- Optional battery charger stand
- Automatic power-off function for power saving



Features



Unique Sensor Design



Bed Rail Mount



Battery Charger Stand



Carrying Bag