

# H10

Finger Oximeter



# H10 Finger Oximeter

H10 Finger Oximeter delivers a fast, reliable oximetry in a very handy, finger-tip solution by combining the monitor and the sensor into one unit. H10 is small in size, nonetheless it provides SpO<sub>2</sub>, pulse rate measurement, bar graph and plethysmogram display.

- Compact & lightweight design
- Dual color OLED display
- SpO<sub>2</sub>, PR, Bar Graph and Plethysmogram display
- One-button operation
- 4-direction display with 6 display modes
- 10-level adjustable brightness
- Low power consumption
- Automatic power off design
- Low-battery indicator
- Standard 2 AAA batteries power supply



## Features



6 Display Modes



Pocket-size Design



One-button Operation

## About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- In-Vitro Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



Edan Instruments, Inc. | No.15 Jinhui Rd., Jinsha Community,  
Kengzi Subdistrict, Pingshan District, Shenzhen | 518122, P.R. China  
+86.755.26898326 | [www.edan.com.cn](http://www.edan.com.cn) | [info@edan.com.cn](mailto:info@edan.com.cn)

© 2018 Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.



ENG-PM-H10-V2.4-20180125